




SOPORTE PARA MARCO CARDIACO PARA CARRO DE PARO



INVENTARIO TÉCNICO

Nombre: SOPORTE PARA MARCO CARDÍACO PARA CARRO DE PARO


Modelo: S/M

Marca: ARMSTRONG MEDICAL

Referencia: ABK-1

Cantidad	Rack	Columna	Tramo
63	E	1	1

VSTRUCTION SHEET ABK-1 CARDIAC BOARD



- 1. Remove the protective covering from the 3 brackets.
- 2. Assemble the brackets to the back of the cabinet with 2 tapping screws each.
- 3. Assemble the rubber bumper with the remaining tapping screw at the left end as shown.

Armstrong Medical Industries Inc.
575 Knightsbridge Parkway • Post Office Box 700 • Lincolnshire, IL 60069-0700
800/323-4220 • Local: 847/913-0101 • Fax: 847/913-0138
www.armstrongmedical.com

BYAM 130